

Kodak

DirectView

SPT1 Radiographic Table



Simplify patient positioning with your KODAK DIRECTVIEW DR System

Positioning patients for digital radiography examinations is convenient, efficient and precise with the mobile KODAK DIRECTVIEW SPT1 Radiographic Table. Designed for maximum manoeuvrability with KODAK DIRECTVIEW DR Systems, this versatile imaging table is also suitable for CR and film-based procedures.

Easy to handle

- ▶ Directional wheel castors and a large handgrip allow for easy movement in any direction.
- ▶ Hydraulic control elevates the table—with no external power source or clumsy cable.
- ▶ This table is ideal for use with the DIRECTVIEW DR 7500 System wall stand.

Floating tabletop

- ▶ A four-way floating tabletop moves easily for precise patient positioning.
- ▶ Foot-pedal control locks the tabletop in position.
- ▶ Tabletop width enhances positioning ease with the DIRECTVIEW DR 3500 System.

Four-wheel brakes

- ▶ A brake on every wheel helps to ensure stability and patient safety.
- ▶ All brakes are set and released simultaneously with one press of a foot pedal. Three pairs of pedals are conveniently located.
- ▶ The table can be anchored to the floor for increased collision avoidance with the DIRECTVIEW DR 9500 System.

Easy to clean

- ▶ Streamlined table design simplifies cleaning and disinfecting.

Maintenance-free

- ▶ This rugged, self-contained table requires no maintenance.

Specifications

Table dimensions

- ▶ Length (including handle): 85.8 in. (218 cm)
- ▶ Width: 28.7 in. (73 cm)
- ▶ Maximum height: 33.9 in. (86 cm)

Tabletop dimensions

- ▶ Length: 83.5 in. (212 cm)
- ▶ Width: 25.6 in. (65 cm)
- ▶ Thickness: 2.1 in. (5.4 cm)

Threshold height

- ▶ Maximum height: 0.6 in. (15 mm)

Table weight

- ▶ 341.7 lbs (155 kg)

Allowable patient weight

- ▶ Maximum allowable patient weight: 500 lbs (225 kg)
- ▶ Lighter than average patient weight: the head of the table will be elevated 20 mm or less

Maximum deflection: 0.79 in. to -1.57 in.
(20 mm to -40 mm)

X-ray transparent surface

- ▶ Length: 83.5 in. (212 cm)
- ▶ Width: 25.6 in. (65 cm)

X-ray absorption factor

- ▶ Factor: ≤ 0.8 mm aluminium equivalence

Travel position of the tabletop

- ▶ Elevating range: 9.8 in. (25 cm)
- ▶ Longitudinal float: 22.5 in. (56 cm)
- ▶ Transversal float: ± 4.8 in. (± 12 cm)

Environmental requirements

Transit and storage

- ▶ The receiving and storage area(s) must be dry and able to provide the proper humidity and temperature control required for the equipment.

Storage temperature

- ▶ -9 to 150°F (-23 to 66°C)

Relative humidity

- ▶ 35 to 75% (non-condensing)

Atmospheric pressure

- ▶ 70 to 106 kPa

Operating temperature

- ▶ 50 to 104°F (10 to 40°C)

Kodak
DirectView
SPT1 Radiographic Table



More information

Learn more about the KODAK DIRECTVIEW SPT1 Radiographic Table and Carestream Health's other products for enhanced DR performance. Contact your Carestream Health representative, or call 1-877-865-6325, ext. 227.

→ carestreamhealth.com

Kodak

Licensed Product

Carestream Health, Inc.
150 Verona Street
Rochester, NY 14608

Carestream 
HEALTH

The innovation in **Kodak** health products